

The AffloVest is reimbursed by Medicare, Medicaid and private insurance for qualified patients.



Patents Pending – USA & Worldwide

AffloVest requires a doctor's prescription for treatment by High Frequency Chest Wall Oscillation (HFCWO). The AffloVest has received the FDA's 510k clearance for U.S. market availability, and is approved for Medicare and private health insurance reimbursement under the Healthcare Common Procedure Coding System (HCPCS) code E0483 – High Frequency Chest Wall Oscillation.

Patients must qualify to meet eligibility requirements.

The AffloVest is the first truly portable, battery operated, fully mobile during use, High Frequency Chest Wall Oscillation (HFCWO) vest. The AffloVest promotes airway clearance and lung secretion mobilization in physician prescribed treatment of respiratory diseases like Bronchiectasis, Cystic Fibrosis, MS, MD, ALS, and other neuromuscular diseases.



Ease of Use

The AffloVest puts complete control in the patient's hands.

- **Fully mobile during use**
- **Digital, programmable controller**
- **Nine total setting variations** let doctors and patients determine the best level of personal treatment for each user
- **Three modes of oscillation treatment**
- **Compliance monitoring**
- **Three adjustable intensity levels** (Soft, Medium, Intense)
- **Quiet during operation**
- **AC/DC power** option allows users to recharge the AffloVest's battery or operate the AffloVest at home or in the car*

* The AffloVest can also be powered using supplied AC adapter plugged into an electrical outlet



No Carts



No loud heavy generators or being tethered to hoses and electrical outlets



"Mobility during use is a huge differentiating factor with this oscillation vest."

– Bill Clark,
Leading COPD Advocate



Freedom of movement utilizing the latest technology

Afflovest is different. It does not use air generator and air bladder technology first developed in the 1980's. No loud and heavy air generators. No air hoses and electrical cords to connect and slow you down. With AffloVest's 21st century technology it gives your doctor a new tool utilizing the latest in oscillation therapy, portability and mobility during use. The AffloVest features built-in oscillating modules that allow for freedom of movement to travel, go for a walk in the park, move around the house, or relax in front of the television. Patient treatment plans can be customized using the targeted motors for specific lung lobes or chest wall zones using different levels of intensity.

- **Rechargeable battery pack** that can last up to 3+ hours on highest settings
- **Afflovest is easy to store and easy to carry;** the Go Anywhere Travel Case makes it perfect for the respiratory patient who travels
- **Available in seven sizes** (XXS, XS, S, M, L, XL, XXL); fits Chest circumference of 18" to 65"+



Targeted treatment, increased patient use, improved lung function.

Doctors recommend and prescribe HFCWO vests for their effectiveness in providing airway clearance. The AffloVest has eight oscillating motors creating eight individual pressure wave forms, providing disruption in the lungs and mobilizing secretions, easing the body's task of discharging harmful excess secretions. These individual motors allow AffloVest to target different areas of the lungs.

The AffloVest offers three distinctly different modes of treatment that can each be set at one of three intensities, for a total of nine possible variations of oscillation treatment. The ability to tailor treatment to a patient's particular needs coupled with ease of use contributes greatly to the AffloVest's effectiveness and increased frequency of use. This can lead to improved lung function as shown in multiple clinician papers.*



7 Sizes Available	Freedom of Movement	9 Setting Variations	Quiet During Operation	Programmable Digital Controller with Compliance Meter	Fully Mobile During Use

Breathe Better. Be Free.
Clear your Chest with AffloVest.™



Patents Pending – USA & Worldwide

INTERNATIONAL
BIOPHYSICS
CORPORATION

2101 E. St. Elmo Rd. Ste. 275 Austin, TX 78744

(T) 888-711-1145 | (F) 888-793-2319 | info@biophysicscorp.com | afflovest.com

MKT0016 Rev. E

The AffloVest is a physician prescribed oscillation treatment device. Respiratory patients should consult their physician to determine how AffloVest can help.

AffloVest[®]
by International Biophysics

**Breathe Better.
Be Free.**



Another innovative product brought to you by

INTERNATIONAL
BIOPHYSICS
CORPORATION

* Cooper, Michael, RT. An evidence-based study of adolescents with cystic fibrosis demonstrated that AffloVest By International Biophysics contributed to improved lung function scores. Chicago, IL.

Tackett, Michelle, RRT, Henderson, Vivian, RRT. Lung function improvement with AffloVest HFCWO use: a clinician's perspective on PFT score data from 25 patients with cystic fibrosis. Knoxville, TN.